

Development of an *in vitro* potency assay for AVB-101, a recombinant AAV9 gene therapy for the treatment of frontotemporal dementia with progranulin mutations



R. Lyth,¹ G. Kooner,¹ M. Donnelly,¹ J. Kenth,¹ I. Blount,¹ A. Bloom¹

AviadoBio Ltd, London, UK

OBJECTIVE

To develop an *in vitro* potency assay which allows for accurate quantification of the progranulin (PGRN) protein to improve and replace an existing *in vivo/in vitro* qualitative potency assay.

CRISPRa – sgRNA guide screen

- Finally, the double-KO, clonal cell line (HEK293 Δ X1 Δ GRN) was further modified to stably express a deactivated Cas9 (dCas9) fused to the hybrid VP64-p65-Rta tripartite activator (VPR), as well as a guide RNA (gRNA) targeting the promoter Y sequence (Figure 4). Binding of the gRNA to the latter sequence elicited transcriptional activation, and thus *hGRN* gene expression from AVB-101.
- Four gRNA sequences targeting promoter Y were screened to determine the single gRNA that brought about the greatest degree of transcriptional activation.
- Guide performance was determined by expression of endogenous *gene Z*, known to be downstream of promoter Y.
- All four gRNA sequences showed an increase in expression of *gene Z*, compared to cells transduced with the negative, non-targeting control (NTC) (Figure 5).
- Single-guide RNA (sgRNA) 2 showed the greatest increase in transcriptional activation, with a 484-fold increase in *gene Z* expression. These cells were selected to develop the AVB-101 potency assay with and referred to as HEK293 Δ X1 Δ GRN-sgRNA2.

Figure 4: Generation of dCAS9-vpr-sgRNA lentiviral stable cell lines





Figure 5: Fold change in expression of *gene Z* normalized to a housekeeping gene (*RPP30*) and non-targeting control cells



RESULTS

- HEK293 Δ X1 Δ GRN-sgRNA2 demonstrated a dose-dependent response in relation to multiplicity of infection (MOI) after transduction with AVB-101 (Figures 6–7).
- No PGRN expression was observed in the NTC cells transduced with AVB-101 (Figure 7).

Figure 7: Transduction of HEK293 Δ X1 Δ GRN-sgRNA2 and NTC cells with AVB-101





Figure 6: Workflow for *in vitro* AVB-101 potency assay




CONCLUSIONS

- We have successfully generated an immortalized engineered cell line with improved AAV9 permissivity, increased sensitivity by reduction of endogenous PGRN expression, and tissue-specific promoter activation.
- The engineered cell line has enabled the development of an *in vitro* potency assay for AVB-101, replacing the previous *in vivo/in vitro* hybrid method.
- The cell line engineering strategy employed provides an elegant solution for gene therapy potency assay development where tissue-specific promoters are utilized.

INTRODUCTION


- AVB-101 (Figure 1) is a recombinant, adeno-associated virus serotype 9 (AAV9)-based gene therapy designed to deliver a functional copy of the human granulin (*hGRN*) gene to cortical neurons, restoring human PGRN to physiological levels and potentially preventing disease progression.^{1,2}
- The AVB-101 therapeutic expression cassette is comprised of the codon optimized *hGRN* driven by a neuronal promoter, denoted as promoter Y and flanked by inverted terminal repeat (ITR) sequences (Figure 1).

Figure 1: AVB-101 key components



- Given the cell type specificity of promoter Y, development of a cell-based potency assay in ubiquitous cell lines, such as HEK293 and HeLa, can be challenging due to negligible expression.
- Whilst neuronal cell lines allow for transcriptional activity, variability arises due to low levels of promoter activity coupled with inconsistent cell growth. Additionally, these cell lines often exhibit a low permissivity to AAV9, resulting in low levels of transgene expression, making dose-dependent detection and quantitation challenging.
- The hybrid *in vivo/in vitro* rat cortical neuron assay (Figure 2) has significant drawbacks; requiring the use of animals and licensed operators, qualitative results, and inconsistencies in cell culture and operator performance. Thus, transfer of this assay to a good manufacturing practice (GMP) environment is difficult.
- Here, we describe the systematic engineering of HEK293 cells using gene editing tools to develop a quantitative *in vitro* potency assay.

Figure 2: Workflow for hybrid *in vivo/in vitro* AVB-101 potency assay




METHODS

Cell line development

- To increase AAV9 transduction, a gene encoding for a critical protein linked to AAV9 binding and cellular uptake (Protein X1) was inactivated by means of CRISPR/Cas9-mediated gene editing (Figure 3A). Several clones were isolated, and a clonal line (HEK293 Δ X1) was generated and assessed for AAV9 transduction (Figure 3B). Δ X1 cells showed improved transduction efficiency (Figure 3B).
- Next, the endogenous *GRN* gene was disrupted in the clonal knockout (KO) cell line (HEK293 Δ X1) to abrogate expression of endogenous PGRN from the cell line, increasing assay sensitivity (Figure 3C). Once again, the modified pool was cloned, and reduced expression of *GRN* in the resulting clone confirmed.

Figure 3: Generation of HEK293 Δ X1 Δ GRN cells



A) Cell line development workflow. B) Peak area of PGRN (identified by automated western blot system [JESS Simple Western™]) is directly related to PGRN concentration. C) In cells where the *GRN* gene was disrupted, endogenous PGRN expression was reduced ~10 fold. Error bars represent standard deviation.

REFERENCES: ¹Lee YB, et al. Oral presentation at ESGCT Annual Congress 2022; ²Miranda CJ, et al. Poster presented at ESGCT Annual Congress 2023.

ABBREVIATIONS: AAV9: adeno-associated virus serotype 9; BSU: biological services and biobanking; Cas9: CRISPR-associated protein 9; CRISPR: clustered regularly interspaced short palindromic repeats; dCas9: deactivated Cas9; ELISA: enzyme-linked immunosorbent assay; GMP: good manufacturing practice; GRN: granulin; gRNA: guide RNA; HEK: human embryonic kidney; hGRN: human granulin; ITR: inverted terminal repeat; KO: knockout; MOI: multiplicity of infection; NTC: non-targeting control; PGRN: progranulin; RT-ddPCR: reverse-transcriptase droplet-digital polymerase chain reaction; sgRNA: single-guide RNA; TTN: titin; VPR: VP64-p65-Rta tripartite activator; WT: wild type.

ACKNOWLEDGMENTS & DISCLOSURES: This study was funded by AviadoBio Ltd. RL, GK, MD, JK, IB and AB are employees and/or shareholders of AviadoBio Ltd. The authors meet criteria for authorship as recommended by the International Committee of Medical Journal Editors (ICMJE). Medical writing and editorial support was provided by Calum Suggett, MSc, of Costello Medical, UK, and funded by AviadoBio Ltd.